



BCAR1 Antibody

Product Code	CSB-PA002593GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56945
Immunogen	Human BCAR1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	breast cancer anti-estrogen resistance 1;BCAR1;CAS;CAS1;CASS1;CRKAS;P130Cas ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BCAR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.