



BMF Antibody

Product Code	CSB-PA002725GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LC9
Immunogen	Human Bmf
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Protein G Purified
Isotype	IgG
Alias	BCL2 modifying factor; Bmf; AW260063; AW411872 ;
Product Type	Purified Rabbit Anti mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BMF
Target Details	<p>This protein belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. This protein contains a single BCL2 homology domain 3 (BH3), and has been shown to bind BCL2 proteins and function as an apoptotic activator. This protein is found to be sequestered to myosin V motors by its association with dynein light chain 2, which may be important for sensing intracellular damage and triggering apoptosis. Alternatively spliced transcript variants encoding different isoforms have been identified.</p>