



BNIP1 Antibody

Product Code	CSB-PA002769GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z465
Immunogen	Human BNIP1
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	BCL2/adenovirus E1B 19kD interacting protein like; BNIP1; BNIP-S; BNIP1-1; BNIP1-2; BNIP1L1; BNIP1L2; PP753 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BNIP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.