



BRCA1 Antibody

Product Code	CSB-PA002794GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38398
Immunogen	Human BRCA1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	breast cancer 1, early onset; BRCA1; BRCAI; BRCC1; BROVCA1; IRIS; PSCP; RNF53 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BRCA1
Target Details	<p>This gene encodes a nuclear phosphoprotein that plays a role in maintaining genomic stability, and it also acts as a tumor suppressor. The encoded protein combines with other tumor suppressors, DNA damage sensors, and signal transducers to form a large multi-subunit protein complex known as the BRCA1-associated genome surveillance complex (BASC). This gene product associates with RNA polymerase II, and through the C-terminal domain, also interacts with histone deacetylase complexes. This protein thus plays a role in transcription, DNA repair of double-stranded breaks, and recombination. Mutations in this gene are responsible for approximately 40% of inherited breast cancers and more than 80% of inherited breast and ovarian cancers. Alternative splicing plays a role in modulating the subcellular localization and physiological function of this gene. Many alternatively spliced transcript variants, some of which are disease-associated mutations, have been described for this gene, but the full-length natures of only some of these variants has been described. A related pseudogene, which is also located on chromosome 17, has been identified.</p>
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.