



AAMDC Antibody

Product Code	CSB-PA003030GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H7C9
Immunogen	Human C11orf67
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	chromosome 11 open reading frame 67;C11orf67;CK067;FLJ21035;MGC3367;PTD015 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AAMDC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.