



C19orf10 Antibody

Product Code	CSB-PA003433GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969H8
Immunogen	Human C19orf10
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	chromosome 19 open reading frame 10;C19orf10;EUROIMAGE1875335;IL25;IL27;IL27w;R33729_1;SF20 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYDGF
Target Details	This protein was previously thought to support proliferation of lymphoid cells and was considered an interleukin. However, this activity has not been reproducible and the function of this protein is currently unknown.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.