



CHTOP Antibody

Product Code	CSB-PA003619GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3Y2
Immunogen	Human C1orf77
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	chromosome 1 open reading frame 77;C1orf77;DKFZp547E1010;MGC131924;MGC86949;RP1-178F15.2;pp7704 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CHTOP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.