



# C4BPA Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA003951GA01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                         |
| <b>Uniprot No.</b>         | P04003  |
| <b>Immunogen</b>           | Human C4BPA   |
| <b>Species Reactivity</b>  | Human   |
| <b>Tested Applications</b> | ELISA, WB, IHC  |
| <b>Form</b>                | Liquid  |
| <b>Storage Buffer</b>      | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles. |
| <b>Purification Method</b> | Antigen Affinity Purified   |
| <b>Isotype</b>             | IgG   |
| <b>Product Type</b>        | Polyclonal Antibody   |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Target Names</b>        | C4BPA   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.           |