



C4BPA Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA003951GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P04003 |
| Immunogen | Human C4BPA |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB, IHC |
| Form | Liquid |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | C4BPA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |