



HAUS3 Antibody

Product Code	CSB-PA003959GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68CZ6
Immunogen	Human C4orf15
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	chromosome 4 open reading frame 15;C4orf15;DKFZp686l1868;IT1;MGC4701 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HAUS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.