



CAB39L Antibody

Product Code	CSB-PA004382GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9S4
Immunogen	Human CAB39L
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	calcium binding protein 39-like;CAB39L;FLJ12577;MO25-BETA;MO2L;RP11-103J18.3;bA103J18.3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CAB39L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.