



CADM1 Antibody

Product Code	CSB-PA004425GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY67
Immunogen	Human CADM1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	cell adhesion molecule 1;CADM1;BL2;DKFZp686F1789;IGSF4;IGSF4A;MGC149785;MGC51880;NEC L2;Necl-2;RA175;ST17;SYNCAM;TSLC1;sTSLC-1;sglGSF;synCAM1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CADM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.