



CALCOCO1 Antibody

Product Code	CSB-PA004436GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P1Z2
Immunogen	Human CALCOCO1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	calcium binding and coiled-coil domain 1;CALCOCO1;Cocoa;KIAA1536;PP13275;calphoglin ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CALCOCO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.