



CAMK2D Antibody

Product Code	CSB-PA004467GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13557
Immunogen	Human CAMK2D
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	calcium/calmodulin-dependent protein kinase II delta;CAMK2D;CAMKD;DKFZp686G23119;DKFZp686I2288;MGC44911;CaMKII delta
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CAMK2D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.