



CAMP Antibody

Product Code	CSB-PA004476GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49913
Immunogen	Human CAMP
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	cathelicidin antimicrobial peptide;CAMP;CAP18;CRAMP;FALL-39;FALL39;HSD26;LL37 ;
Immunogen Species	Homo sapiens (Human)
Target Names	CAMP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.