



CBLL1 Antibody

Product Code	CSB-PA004581GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q75N03
Immunogen	Human CBLL1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Cas-Br-M (murine) ecotropic retroviral transforming sequence-like 1;CBLL1;FLJ23109;HAKAI;MGC163401;MGC163403;RNF188 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CBLL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.