



CD109 Antibody

Product Code	CSB-PA004878GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6YHK3
Immunogen	Human CD109
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	CD109 molecule; CD109; CPAMD7; DKFZp762L1111; FLJ38569; FLJ41966; RP11-525G3.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD109
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.