



CD274 Antibody

Product Code	CSB-PA004911GA01HU
Abbreviation	CD274
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZQ7
Immunogen	Human CD274
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IP, IHC, IF, FC, ELISA
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	CD274 molecule;CD274;B7-H;B7H1;MGC142294;MGC142296;PD-L1;PDCD1L1;PDCD1LG1;PDL1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Species	Homo sapiens (Human)
Target Names	CD274