



CD36 Antibody

Product Code	CSB-PA004927GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16671
Immunogen	Human CD36
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,IHC
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	CD36 molecule (thrombospondin receptor);CD36;CHDS7;FAT;GP3B;GP4;GPIV;PASIV;SCARB3 ;
Immunogen Species	Homo sapiens (Human)
Target Names	CD36
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.