



CEACAM21 Antibody

Product Code	CSB-PA005162GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3KPI0
Immunogen	Human CEACAM21
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	carcinoembryonic antigen-related cell adhesion molecule 21;CEACAM21;CEACAM3;FLJ13540;MGC119874;R29124_1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CEACAM21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.