



CMPK1 Antibody

Product Code	CSB-PA005616GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30085
Immunogen	Human CMPK1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	cytidine monophosphate (UMP-CMP) kinase 1, cytosolic;CMPK1;CMK;CMPK;UMK;UMP-CMPK;UMPK ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CMPK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.