



CPA3 Antibody

Product Code	CSB-PA005877GA01HU
Abbreviation	CPA3
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15088
Immunogen	Human CPA3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	carboxypeptidase A3 (mast cell);
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Species	Homo sapiens (Human)
Target Names	CPA3
Target Details	Three different forms of human pancreatic procarboxypeptidase A have been isolated. This gene encodes a form which is obtained as a binary complex of a procarboxypeptidase A with proproteinase E and functions as a secretory granule metalloexopeptidase.