



CPT1C Antibody

Product Code	CSB-PA005924GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCG5
Immunogen	Human CPT1C
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	carnitine palmitoyltransferase 1C;CPT1C;CATL1;CPT1P;CPTIC;FLJ23809 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CPT1C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.