



# CRISP3 Antibody

<b>Product Code</b>	CSB-PA005976GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P54108
<b>Immunogen</b>	Human CRISP3
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA,WB,IHC
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	cysteine-rich secretory protein 3;CRISP3;Aeg2;CRISP-3;CRS3;MGC126588;SGP28;dJ442L6.3 ;
<b>Product Type</b>	Purified Rabbit Anti Human,Mouse,Rat PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CRISP3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.