



CSNK1A1L Antibody

Product Code	CSB-PA006066GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N752
Immunogen	Human CSNK1A1L
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	casein kinase 1, alpha 1-like;CSNK1A1L;MGC33182;RP11-532O21.2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CSNK1A1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.