



CTAGE1 Antibody

Product Code	CSB-PA006128GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RT6
Immunogen	Human CTAGE1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	cutaneous T-cell lymphoma-associated antigen 1;CTAGE1;CTAGE-1;CTAGE-2 ;CTAGE; CT21.1; CT21.2,CTAGE2
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CTAGE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.