



CTNNA2 Antibody

Product Code	CSB-PA006165GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26232
Immunogen	Human CTNNA2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	catenin (cadherin-associated protein), alpha 2;CTNNA2;CAP-R;CAPR;CTNR;DKFZp686H02198 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CTNNA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.