



CXCL14 Antibody

Product Code	CSB-PA006244GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95715
Immunogen	Human CXCL14
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	chemokine (C-X-C motif) ligand 14;CXCL14;BMAC;BRAK;KS1;Kec;MGC10687;MIP-2g;NJAC;SCYB14;bolekine ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CXCL14
Target Details	This gene belongs to the cytokine gene family which encode secreted proteins involved in immunoregulatory and inflammatory processes. This protein is structurally related to the CXC (Cys-X-Cys) subfamily of cytokines. Members of this subfamily are characterized by two cysteines separated by a single amino acid. This cytokine displays chemotactic activity for monocytes but not for lymphocytes, dendritic cells, neutrophils or macrophages. It has been implicated that this cytokine is involved in the homeostasis of monocyte-derived macrophages rather than in inflammation.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.