



CXCL3 Antibody

Product Code	CSB-PA006249GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19876
Immunogen	Human CXCL3
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	chemokine (C-X-C motif) ligand 3;CXCL3;CINC-2b;GRO3;GROg;MIP-2b;MIP2B;SCYB3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CXCL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.