



DCT Antibody

Product Code	CSB-PA006562GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40126
Immunogen	Human DCT
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	dopachrome tautomerase (dopachrome delta-isomerase, tyrosine-related protein 2);DCT;TYRP2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DCT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.