



DDX51 Antibody

Product Code	CSB-PA006632GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N8A6
Immunogen	Human DDX51
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	DEAD (Asp-Glu-Ala-Asp) box polypeptide 51;DDX51;DKFZp686N2081;MGC42193 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DDX51
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.