



DEF6 Antibody

Product Code	CSB-PA006650GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4E7
Immunogen	Human DEF6
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	differentially expressed in FDCP 6 homolog (mouse);DEF6;IBP;SLAT ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DEF6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.