



DIAPH3 Antibody

Product Code	CSB-PA006892GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NSV4
Immunogen	Human DIAPH3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	diaphanous homolog 3 (Drosophila);DIAPH3;DKFZp434C0931;DKFZp686A13178;DRF3;FLJ34705;dia p3;mDia2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DIAPH3