



DIS3 Antibody

Product Code	CSB-PA006911GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2L1
Immunogen	Human DIS3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	DIS3 mitotic control homolog (S. cerevisiae);DIS3;DKFZp667L1817;EXOSC11;FLJ10484;KIAA1008;MGC33035;RP11-342J4.3;RRP44;bA555G22.1;dis3p ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DIS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.