



DLEU1 Antibody

Product Code	CSB-PA006930GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43261
Immunogen	Human DLEU1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	deleted in lymphocytic leukemia 1 (non-protein coding);DLEU1;BCMS;DLB1;DLEU2;LEU1;LEU2;MGC22430;NCRNA00021;XT P6 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DLEU1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.