



DLGAP5 Antibody

Product Code	CSB-PA006944GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15398
Immunogen	Human DLGAP5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	discs, large (Drosophila) homolog-associated protein 5; DLGAP5; DLG7; HURP; KIAA0008 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DLGAP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.