



DMKN Antibody

Product Code	CSB-PA006965GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6E0U4
Immunogen	Human DMKN
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	dermokine;DMKN;FLJ79374;FLJ79374,;UNQ729;ZD52F10 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DMKN
Target Details	This gene is upregulated in inflammatory diseases, and it was first observed as expressed in the differentiated layers of skin. The most interesting aspect of this gene is the differential use of promoters and terminators to generate isoforms with unique cellular distributions and domain compliments. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene, but the full-length nature of some of them has not been determined.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.