



DNAJB11 Antibody

Product Code	CSB-PA007004GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBS4
Immunogen	Human DNAJB11
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	DnaJ (Hsp40) homolog, subfamily B, member 11;DNAJB11;ABBP-2;ABBP2;DJ9;EDJ;ERdj3;ERj3;HEDJ;PRO1080;UNQ537;hDj9 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DNAJB11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.