



DNAJC10 Antibody

Product Code	CSB-PA007017GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXB1
Immunogen	Human DNAJC10
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	DnaJ (Hsp40) homolog, subfamily C, member 10;DNAJC10;DKFZp434J1813;ERdj5;JPDI;MGC104194 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DNAJC10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.