



DOCK9 Antibody

Product Code	CSB-PA007105GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZ29
Immunogen	Human DOCK9
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	dedicator of cytokinesis 9;DOCK9;DKFZp686C11110;DKFZp686D2047;DKFZp686N04132;FLJ11949;FLJ16744;FLJ44528;FLJ45282;FLJ45601;KIAA1058;KIAA1085;ZIZ1;ZIZIMIN1 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DOCK9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.