



DYNC1I2 Antibody

Product Code	CSB-PA007294GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13409
Immunogen	Human DYNC1I2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	dynein, cytoplasmic 1, intermediate chain 2;DYNC1I2;DNCI2;FLJ21089;FLJ90842;IC2;MGC104199;MGC111094;MGC9324 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DYNC1I2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.