



EDN1 Antibody

Product Code	CSB-PA007400GA01HU
Abbreviation	EDN1
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05305
Immunogen	Human EDN1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	endothelin 1;EDN1;ET1;HDLCQ7 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Species	Homo sapiens (Human)
Target Names	EDN1
Target Details	This protein is proteolytically processed to release a secreted peptide termed endothelin 1. This peptide is a potent vasoconstrictor and is produced by vascular endothelial cells. Endothelin 1 also can affect the central nervous system. Two transcript variants encoding different isoforms have been found for this gene.