



# CRACR2A Antibody

|                     |                                                                                         |
|---------------------|-----------------------------------------------------------------------------------------|
| Product Code        | CSB-PA007444GA01HU                                                                      |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                           |
| Uniprot No.         | Q9BSW2                                                                                  |
| Immunogen           | Human EFCAB4B                                                                           |
| Raised In           | Rabbit                                                                                  |
| Species Reactivity  | Human                                                                                   |
| Tested Applications | ELISA, WB                                                                               |
| Storage Buffer      | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.   |
| Purification Method | Antigen Affinity purified                                                               |
| Isotype             | IgG                                                                                     |
| Alias               | EF-hand calcium binding domain<br>4B;EFCAB4B;DKFZp686G13246;FLJ33046;FLJ33805;MGC4266 ; |
| Product Type        | Purified Rabbit Anti human PolyClonal Antibody                                          |
| Immunogen Species   | Homo sapiens (Human)                                                                    |
| Target Names        | CRACR2A                                                                                 |
| Usage               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.             |