



EID1 Antibody

Product Code	CSB-PA007499GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6B2
Immunogen	Human EID1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	EP300 interacting inhibitor of differentiation 1; EID1; C15orf3; CR11; EID-1; IRO45620; MGC138883; MGC138884; PNAS-22; PT D014; RBP21 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EID1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.