



# EIF1AX Antibody

<b>Product Code</b>	CSB-PA007506GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P47813
<b>Immunogen</b>	Human EIF1AX
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	eukaryotic translation initiation factor 1A, X-linked;EIF1AX;EIF1A;EIF1AP1;EIF4C;eIF-1A;eIF-4C ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	EIF1AX
<b>Target Details</b>	This gene encodes an essential eukaryotic translation initiation factor. The protein is required for the binding of the 43S complex (a 40S subunit, eIF2/GTP/Met-tRNAi and eIF3) to the 5 end of capped RNA.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.