



EIF2AK2 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA007511GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P19525 |
| Immunogen | Human EIF2AK2 |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA,WB,IHC,IF |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | eukaryotic translation initiation factor 2-alpha kinase 2;EIF2AK2;EIF2AK1;MGC126524;PKR;PRKR ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | EIF2AK2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |