



EIF2S1 Antibody

Product Code	CSB-PA007523GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05198
Immunogen	Human EIF2S1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	eukaryotic translation initiation factor 2, subunit 1 alpha, 35kDa;EIF2S1;EIF-2;EIF-2A;EIF-2alpha;EIF2;EIF2A ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF2S1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.