



EIF3H Antibody

Product Code	CSB-PA007537GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15372
Immunogen	Human EIF3H
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	eukaryotic translation initiation factor 3, subunit H;EIF3H;EIF3S3;MGC102958;eIF3-gamma;eIF3-p40 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF3H
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.