



EIF3J Antibody

Product Code	CSB-PA007539GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75822
Immunogen	Human EIF3J
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	eukaryotic translation initiation factor 3, subunit J;EIF3J;EIF3S1;eIF3-alpha;eIF3-p35 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF3J
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.