



ELF2 Antibody

Product Code	CSB-PA007597GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15723
Immunogen	Human ELF2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	E74-like factor 2 (ets domain transcription factor);ELF2;EU32;NERF;NERF-1A;NERF-1B;NERF-1a,b;NERF-2 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ELF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.